



PROTEGO[®]

RADIATION PROTECTION SYSTEM





PROTEGO[®]

Comprehensive Shielding

BOXES OUT

Harmful Radiation Exposure



>99%
reduction in RE*

*based on the most recent clinical study
from Cardiovascular Revascularization Medicine

all studies are available to read on our website

 www.getprotego.com

How the Protego System Protects You

Protego provides greater than 99% radiation reduction based on real-time dosimetry

Best in class, the Protego Radiation Protection System offers the protection you need to go Apron-Free™*.

Protego provides the largest area of protection


Other systems protect the table area, while the Protego System offers a broad zone of protection as illustrated in the radiation footprint mapped by a facility's physicist.

Protego pans with your table and fits in your current workflow

The Protego System is passive and becomes a part of your table.

Protego's flexible design allows for steep angles in caudal, cranial, and LAO/RAO

The Protego Shield's shape and flexible attachment allow for angles 45° and greater simply by moving your C-arm into position.



idi
imagediagnostics.com

INTELLIGENT DESIGN ► TANGIBLE OUTCOMES



Imagine APRON-FREE IMAGING for Your Whole Team!



Immediate emergent removal

When you need immediate access to the patient, a quick tug upward releases the Protego Shield to easily move it out of your way. See a video demonstration on our website.

Here's what users are saying:

"I've taken the equivalence of 30,000 chest x-rays to my face."

"I don't want any more."

David G. Rizik MD

Radial, jugular, and groin access

The Protego Shield provides easy access to your patient.

An open operating area, free from obstruction

Nothing comes between you and your access site, allowing for better visualization and posture while you work. The Protego System stays out of your way for better table side ergonomics.

Low cost per procedure

With re-usable scatter shields and only three drapes we keep your per case cost low.

Set-up less than 3 minutes, using only 3 drapes

Most of the system stays on the table for fast and easy set up.

Excellent patient visualization and communication

The healthcare team can keep a close and clear eye on the patient at all times with our integral dual-camera and monitor system.

All with a small footprint

Most of the system is comprised of shields residing on the table, the above-table Protego Shield, and the Mobile Shield to extend radiation protection to greater area. The Protego Shield stand uses the same floor space as a monitor cart.

"The biggest question for everyone in the lab was how much radiation exposure would be reduced with Protego. My RaySafe radiation badge registered zero.

It was the same for all staff badges on my side of the shielding wall."

"We performed "lead-free cases" all week and the near-zero radiation reading in the shielded zone were impressively low: near zero in every case."

Mort Kern, MD

Procedure parameters for Apron-Free™ imaging at your facility must be defined by your Radiation Safety Officer or facility physicist as well as comply with state and local regulations.

INTELLIGENT DESIGN ► TANGIBLE OUTCOMES



ORDER: sales@imagediagnostics.com
SUPPORT: techsupport@imagediagnostics.com
WEBSITE: <https://imagediagnostics.com>

+1-978-829-0009
Monday through Friday
8am to 5pm, Eastern

